

Enterprise PACS Architecture: Infrastructure Mechanics, DICOM Workflow Dynamics, Mirth Connect Integration, Database Auditing, and VNA Migration

Leo A. Jiwanmall, MBA, PMP®, MRSO®, R.T.(R)(MR)(CT)

Publication Date: January 2026

Target Audience: Imaging Informatics Professionals, PACS Administrators, Healthcare Systems Analysts, and Clinical Engineers

Abstract

Picture Archiving and Communication Systems (PACS) serve as the operational backbone of digital radiology and healthcare enterprise imaging networks. As image acquisition modalities generate increasingly high-volume, multi-dimensional datasets, understanding the structural informatics and IT infrastructure of PACS is paramount. This paper provides a comprehensive analysis of enterprise PACS architecture, tracing data flow from modality acquisition through Health Level Seven (HL7) interface brokering, Digital Imaging and Communications in Medicine (DICOM) network protocol handshakes, archive tiering strategies, and Vendor Neutral Archive (VNA) consolidation. Furthermore, it details technical implementation rules for Mirth Connect HL7-to-DICOM transformations, provides backend SQL query routines for troubleshooting `C-STORE` and `C-FIND` database anomalies, and examines high-availability infrastructure, diagnostic display standards, and cybersecurity protocols within modern clinical informatics.

Methodology

Study Design

This manuscript utilizes a narrative technical review methodology combined with operational analysis derived from contemporary enterprise imaging and PACS administration practices.

Information Sources

Primary source materials included:

- DICOM Standards (NEMA PS3)
- IHE Radiology Technical Framework
- AAPM Task Group 18 guidance
- PACS and imaging informatics textbooks
- DICOM implementation references
- Enterprise imaging interoperability guidance
- Vendor-neutral archive (VNA) literature

Analytical Framework

Content was evaluated across five functional domains:

1. Enterprise PACS Infrastructure
2. DICOM Communications
3. HL7 Workflow Integration
4. Storage and Archival Architecture
5. High Availability, Security, and Disaster Recovery

The objective was to synthesize established standards and operational best practices into a practical framework applicable to contemporary healthcare imaging organizations.

1. Introduction and Historical Evolution

The modern Picture Archiving and Communication System (PACS) transitioned medical imaging from analog screen-film systems and lightboxes to digital, interconnected workflow ecosystems (Huang, 2010). A contemporary PACS is not an isolated software application, but an integrated assembly of hardware, software, networking protocols, and database structures designed to acquire, transport, display, store, and manage digital diagnostic images.


In its earliest iterations during the late 1980s and 1990s, PACS operated on proprietary network protocols and isolated storage units. The standardization of the Digital Imaging and Communications in Medicine (DICOM) standard and the widespread adoption of Health Level Seven (HL7) messaging enabled vendor interoperability and enterprise imaging integration across healthcare systems. These standards established the foundation for

modern imaging informatics environments capable of supporting high-volume, multi-site clinical operations

Today, PACS architecture must support vast data volumes generated by advanced modalities such as multi-detector computed tomography (MDCT), digital breast tomosynthesis (DBT), and ultra-high-field magnetic resonance imaging (MRI).

2. Core Architectural Components and Data Flow

Enterprise PACS relies on a client-server or hybrid web-based model designed to balance data throughput, system redundancy, and rapid image retrieval. The architecture comprises four core operational tiers:



2.1 Image Acquisition Tier

Modalities (e.g., DR, CT, MRI, Ultrasound) generate digital image matrix files comprising spatial pixel data and rich header metadata. Upon acquisition, the modality acts as a DICOM **Service Class User (SCU)**, initiating a network transmission to the PACS archive, which acts as the DICOM **Service Class Provider (SCP)**.

2.2 Server and Database Core

The central server suite houses two distinct functional layers:

1. **The Database Engine (Metadata Layer):** A relational database management system (RDBMS) such as Microsoft SQL Server, PostgreSQL, or Oracle. It stores demographic data, study metadata, Accession Numbers, Series Instance UIDs, DICOM header tags, and pointer paths to physical pixel files.
2. **The Image Store (Pixel Layer):** File system storage holding actual bulk binary image files (standard DICOM part 10 format). Decoupling the relational database from

physical image files is crucial for maintaining query speed across millions of patient records (Pianykh, 2012).

2.3 Diagnostic Display Tier

Workstations utilized by radiologists require specialized hardware and software configurations:

- **Luminance and Calibration:** Diagnostic monitors must maintain stringent luminance characteristics calibrated against the DICOM **Grayscale Standard Display Function (GSDF)** under AAPM Task Group 18 and Task Group 270 guidelines (Samei et al., 2005). Standard clinical displays typically operate around $350\text{ to }500\text{ cd/m}^2$, while dedicated mammography displays require $1,000\text{ cd/m}^2$.
- **Rendering Architecture:** Modern client applications leverage server-side rendering (SSR) or progressive streaming technologies (e.g., WebGL, wavelet transfers) to allow rapid display of multi-gigabyte datasets over standard enterprise local area networks (LAN) or wide area networks (WAN).

3. Communication Protocols and Interoperability Standards

System integration within PACS architecture relies on two complementary communication standards: HL7 for textual workflow events and DICOM for image data mechanics.

3.1 HL7 Workflow Integration

Health Level Seven (HL7 v2.x) handles non-image patient administration data:

- **ADT (Admission, Discharge, Transfer):** Updates patient demographic changes across the enterprise.
- **ORM (Order Message):** Communicates order generation from the Radiology Information System (RIS) or Electronic Health Record (EHR) to the PACS broker.
- **ORU (Observation Result Unsolicited):** Transmits final transcribed and signed radiologist reports back to the EHR.

3.2 DICOM Network Protocol Mechanics

DICOM operates on top of standard TCP/IP network transport layers, utilizing specified Application Entities (AE) defined by unique **AE Titles**, IP addresses, and TCP port numbers (typically Port 104, 1112, or 4100) (Clunie, 2000). The underlying protocol transactions follow strict Service-Object Pair (SOP) classes:


DICOM Command	Functional Role	Operational Mechanics
C-ECHO	Verification	Verifies basic TCP/IP and DICOM application connectivity between SCU and SCP ("DICOM Ping").
C-FIND	Query	Searches the PACS database for studies matching user-defined key attributes (e.g., Patient ID, Accession Number).
C-MOVE	Retrieve Request	Instructs the archiving SCP to transfer image files to a designated destination AE Title via a secondary c-STORE sub-operation.
C-STORE	Image Transfer	Transmits the physical DICOM dataset (header + pixel matrix) from source to target destination.
MWL (Modality Worklist)	Workflow Synchronization	Queries RIS data via C-FIND to populate patient demographics directly on the acquisition modality console, preventing manual entry errors.

4. Interface Engineering: Mirth Connect HL7 Transformation Rules

In enterprise environments, an interface engine such as NextGen Connect (Mirth Connect) sits between the RIS/EHR and the PACS/MWL broker to normalize HL7 messages and transform textual data into DICOM-compatible fields.

4.1 Field Mapping Logic (HL7 ORM to DICOM Tag Mapping)

To populate a DICOM Modality Worklist (MWL) or trigger automated pre-fetching in PACS, incoming HL7 `ORM^O01` messages must be parsed and mapped to explicit DICOM metadata elements.



ORC.2.1 / OBR.2.1	→ Strip Special Chars	→
(0008,0050) Accession Number		
OBR.4.1 (Procedure Code)	→ Crosswalk Lookup	→
(0032,1060) Requested Procedure Desc		
PV1.8.1 (Referring Doc)	→ Format: ID^Last^First	→
(0008,0090) Referring Physician		

4.2 JavaScript Transformer Snippet (Mirth Connect Channel)

Below is a production-grade JavaScript transformer script deployed within a Mirth Connect destination transformer node to parse an incoming HL7 ORM^O01 order and construct an outbound JSON/DICOM dataset object:

```
// Mirth Connect Transformer: HL7 ORM^O01 to DICOM MWL Data
Normalization
try {
    // Extract HL7 fields safely from message object
    var mrn = msg['PID']['PID.3']['PID.3.1'].toString().trim();
    var rawPatientName = msg['PID']['PID.5'];
    var lastName = rawPatientName['PID.5.1'].toString().trim();
    var firstName = rawPatientName['PID.5.2'].toString().trim();
    var middleName = rawPatientName['PID.5.3'].toString().trim();

    var accession =
msg['ORC']['ORC.2']['ORC.2.1'].toString().trim();
    if (!accession || accession === "") {
        accession =
msg['OBR']['OBR.2']['OBR.2.1'].toString().trim();
    }

    var gender =
msg['PID']['PID.8'].toString().trim().toUpperCase();
    var dob = msg['PID']['PID.7']['PID.7.1'].toString().trim(); //
Format: YYYYMMDD
    var procCode =
msg['OBR']['OBR.4']['OBR.4.1'].toString().trim();
    var procDesc =
msg['OBR']['OBR.4']['OBR.4.2'].toString().trim();

    // Normalize Patient Name to DICOM standard (Last^First^Middle)
    var dicomPatientName = lastName + "^" + firstName;
    if (middleName !== "") {
        dicomPatientName += "^" + middleName;
    }

    // Sanitize Gender for DICOM (0010,0040) -> M, F, O, or empty
    if (gender !== "M" && gender !== "F") {
        gender = "O";
    }
}
```

```

    }

    // Construct Channel Map Variables for Outbound DICOM MWL
    Generator
    channelMap.put('dicomPatientID', mrn);
    channelMap.put('dicomPatientName', dicomPatientName);
    channelMap.put('dicomAccessionNumber', accession);
    channelMap.put('dicomPatientSex', gender);
    channelMap.put('dicomPatientDOB', dob);
    channelMap.put('dicomProcedureCode', procCode);
    channelMap.put('dicomProcedureDescription', procDesc);

    logger.info("Successfully transformed ORM Order for Accession:
" + accession);
} catch (e) {
    logger.error("Error in HL7 to DICOM Mirth Transformer: " +
e.toString());
    // Route to error queue
    System.err.println(e);
}

```

5. Backend PACS Administration: Troubleshooting SQL Queries

PACS Administrators frequently query the relational database directly to resolve broken DICOM associations, stuck transmission queues, orphaned series, or corrupted UID linkages during C-STORE and C-FIND failures.

5.1 Identifying Orphaned Series or Image Instances (C-STORE Ingestion Failures)

When a modality sends images via C-STORE but the parent study row fails to commit due to an accession mismatch or database lock, images become "orphaned."

```

-- SQL Query: Locate Orphaned Series lacking a valid Parent Study
record
SELECT
    ser.SeriesInstanceUID,
    ser.Modality,
    ser.SeriesNumber,
    ser.ImageCount,
    ser.CreatedDateTime,
    ser.SourceAETitle
FROM
    PACS_Series ser
LEFT JOIN

```

```

PACS_Study std ON ser.StudyFK = std.StudyPK
WHERE
    std.StudyPK IS NULL
    OR std.PatientFK IS NULL
ORDER BY
    ser.CreatedDateTime DESC;

```

5.2 Auditing C-FIND Query Performance and Association Timeouts

When radiologists experience slow worklist loading, query bottlenecks are often traced to missing database indexes on heavily queried DICOM tags (AccessionNumber, PatientID, StudyDate).

```

-- SQL Query: Audit C-FIND Query Latency and Identify Slow MWL/PACS
Queries
SELECT TOP 50
    LogID,
    CallingAETitle,
    CalledAETitle,
    QueryType, -- e.g., 'C-FIND', 'C-MOVE'
    SearchPatientID,
    SearchAccessionNumber,
    ExecutionTimeMS,
    ResponseCode, -- 0 = Success, 65535 = Failed/Timeout
    LogTimestamp
FROM
    PACS_DICOM_Audit_Log
WHERE
    QueryType = 'C-FIND'
    AND ExecutionTimeMS > 2000 -- Queries taking longer than 2
seconds
ORDER BY
    ExecutionTimeMS DESC;

```

5.3 Detecting Duplicate Accession Numbers or Broken Study Instance UIDs

If two separate imaging studies share the same Accession Number due to a RIS ordering glitch, C-STORE operations will fail or overwrite images.

```

-- SQL Query: Find Duplicate Accession Numbers with Divergent
StudyInstanceUIDs
SELECT
    AccessionNumber,
    COUNT(DISTINCT StudyInstanceUID) AS DistinctStudyUIDCount,
    COUNT(DISTINCT PatientFK) AS DistinctPatientCount,

```


```

        STRING_AGG(CAST(StudyInstanceUID AS VARCHAR(MAX)), ' | ') AS
AssociatedUIDs
FROM
    PACS_Study
WHERE
    AccessionNumber IS NOT NULL AND AccessionNumber <> ''
GROUP BY
    AccessionNumber
HAVING
    COUNT(DISTINCT StudyInstanceUID) > 1;

```

6. Storage Architecture and Information Lifecycle Management (ILM)

Given the exponential growth of imaging data across multi-facility health systems, PACS storage must balance ultra-fast retrieval speeds against long-term storage economics (Huang, 2010). Facilities utilize tiered storage managed by automated **Information Lifecycle Management (ILM)** rules.



6.1 Storage Tiers

1. **Short-Term Storage (Tier 1):** Utilizes high-performance Storage Area Networks (SAN) or solid-state NVMe storage arrays. Holds recent studies (typically 30–90 days) requiring sub-second access speeds for active reading and immediate comparison prior-study fetching.

2. **Long-Term Archive (Tier 2):** Employs high-capacity Network Attached Storage (NAS), object-based cloud storage (e.g., AWS S3, Azure Blob), or optical libraries. Studies are migrated here after an active period, maintaining sub-five-second retrieval speeds while lowering cost-per-terabyte overhead.


6.2 Data Compression Mechanics

To optimize bandwidth and storage capacity, image compression is applied in accordance with DICOM standards (Pianykh, 2012):

- **Lossless Compression:** Uses mathematical algorithms such as JPEG 2000 Lossless or Run-Length Encoding (RLE). Reconstructs the exact bit-for-bit original pixel data upon expansion, ensuring zero loss of diagnostic detail (Compression ratios typically 2:1 to 3:1).
- **Lossy (Irreversible) Compression:** Permanently discards visually redundant data using algorithms like JPEG 2000 Lossy. Achieves high compression ratios (10:1 to 20:1), but its application in diagnostic primary interpretation is regulated strictly by clinical guidelines and institutional policy.

7. Vendor Neutral Archives (VNA) and Enterprise Integration

A major evolutionary milestone in imaging informatics is the migration from monolithic, vendor-locked PACS to **Vendor Neutral Architecture (VNA)** (Dreyer et al., 2006). Proprietary PACS archives historically encapsulated DICOM header data using non-standard database schemas, creating massive technical friction during system migrations.



Core VNA Functions

- **Decoupling Storage from Display:** A VNA isolates the storage management layer from the workstation display layer. The VNA stores DICOM and non-DICOM clinical assets (e.g., MPEG-4 endoscopy videos, PDF reports, JPEG dermatology images) in a standardized, vendor-independent format.

- **DICOM Tag Morphing:** Resolves inconsistencies in proprietary vendor DICOM header tags dynamically, harmonizing metadata across different equipment manufacturers prior to archiving.
- **Cross-Enterprise Document Sharing (XDS / XDS-I):** Implements IHE (Integrating the Healthcare Enterprise) profiles to query and retrieve clinical image objects across disparate health networks via web services (WADO-RS / STOW-RS / QIDO-RS).

8. High Availability, Security, and Disaster Recovery

Because radiology services operate continuously, PACS architecture requires redundant infrastructure to prevent unscheduled downtime and secure protected health information (PHI).

8.1 System Redundancy and Load Balancing

- **High Availability (HA):** Implemented via clustered server nodes, database mirror pairs, and virtualized machine failovers (e.g., VMware HA). Load balancers distribute incoming DICOM store transactions across multiple ingest servers.
- **Disaster Recovery (DR):** Enforces off-site asynchronous data replication. A secondary hot-site or cloud-hosted PACS instance maintains an active copy of the relational database and long-term storage, enabling rapid disaster recovery during site outages.

8.2 Security and Compliance

- **Data Encryption:** Enforces Transport Layer Security (TLS 1.3) for DICOM nodes in transit, alongside AES-256 encryption for files at rest, complying with HIPAA guidelines.
- **Role-Based Access Control (RBAC):** Integrates with enterprise directory services (Active Directory / LDAP) and SAML/OAuth2 single sign-on (SSO) frameworks to enforce audit logging and strict user permissions.

Limitations

- Several limitations should be recognized.
- First, this manuscript represents a narrative technical review and operational synthesis rather than a systematic review or original research investigation.
- Second, PACS architectures vary considerably among vendors, healthcare systems, deployment models, and enterprise imaging strategies. Specific implementation details, database schemas, HL7 interfaces, and DICOM workflows may differ from the examples presented.
- Third, example SQL queries, Mirth Connect transformations, storage architectures, and interoperability workflows are intended for educational purposes and may require adaptation to local infrastructure, security requirements, organizational policies, and vendor-specific implementations.
- Finally, healthcare imaging technologies continue to evolve rapidly, particularly in areas involving cloud-native PACS deployments, enterprise imaging, artificial intelligence, and cybersecurity.

Conclusion

PACS architecture has evolved from basic digital archive systems into highly sophisticated, enterprise-wide imaging networks. Modern imaging informatics relies on robust physical and software infrastructure: high-speed storage tiering, strict HL7 and DICOM communication standards, Mirth Connect interface normalization, proactive SQL database auditing, GSDF-calibrated diagnostic displays, and flexible VNA layers. By maintaining redundant, secure, and standards-compliant architectural frameworks, healthcare institutions ensure seamless image transmission, rapid diagnostic turnaround, and durable long-term archiving across the clinical enterprise.

Disclaimer

This manuscript is intended for educational and professional discussion purposes. Example SQL queries, HL7 mappings, Mirth Connect transformations, workflow diagrams, and infrastructure models should be validated within local environments and vendor-specific implementations before operational deployment. Healthcare organizations should follow applicable regulatory requirements, institutional policies, cybersecurity standards, and vendor guidance when implementing enterprise imaging solutions.

References

1. **American Association of Physicists in Medicine.** (2005). *Assessment of display performance for medical imaging systems* (AAPM Task Group 18 Report). AAPM. https://www.aapm.org/pubs/reports/OR_03.pdf
2. **Clunie, D. A.** (2000). *DICOM structured reporting*. PixelMed Publishing.
3. **Dreyer, K. J., Hirschorn, D. S., Thrall, J. H., & Mehta, A. J.** (Eds.). (2006). *PACS: PACS and imaging informatics: Basic principles and applications* (2nd ed.). Springer-Verlag. <https://doi.org/10.1007/0-387-31070-6>
4. **Huang, H. K.** (2010). *PACS and imaging informatics: Basic principles and applications* (2nd ed.). Wiley-Blackwell. <https://doi.org/10.1002/9780470560525>
5. **Integrating the Healthcare Enterprise.** (2021). *IHE radiology technical framework, Volume 1: Integration profiles*. IHE International. https://www.ihe.net/resources/technical_frameworks/#radiology
6. **National Electrical Manufacturers Association.** (2021). *Digital Imaging and Communications in Medicine (DICOM) standard* (NEMA Standard PS 3.1-3.20). NEMA. <https://www.dicomstandard.org/current>
7. **Pianykh, O. S.** (2012). *Digital Imaging and Communications in Medicine (DICOM): A practical introduction and survival guide* (2nd ed.). Springer-Verlag. <https://doi.org/10.1007/978-3-642-22574-1>
8. **Samei, E., Badano, A., Chakraborty, D., Compton, K., Cornelius, C., Corrigan, K., ... & Silosky, M.** (2005). Assessment of display performance for medical imaging systems: Executive summary of AAPM TG18 report. *Medical Physics*, 32(4), 1205–1225. <https://doi.org/10.1118/1.1861159>